


FLAIR SLT™ Buccal Tube

GREAT....JUST GOT BETTER!

Introducing the FLAIR SLT™ Buccal Tube, the same robust self-ligating clip of the FLAIR SLT™ Bracket is now available in a buccal tube. Designed specifically to withstand the rigors of treatment and continues to offer easy, uncomplicated archwire changes now through to the molars.



Roth prescription

Self-ligating Buccal Tube	Torque	Ang	Distal Offset	U - R .018	U - L .018	U - R .022	U - L .022
1. & 2. Molar - bondable	-10°	0°	0°	205M-16/27	205M-16/27	255M-16/27	255M-16/27
1. & 2. Molar extended base - bondable	-10°	0°	0°	205M-16/27E	205M-16/27E	255M-16/27E	255M-16/27E
Self-ligating Buccal Tube	Torque	Ang	Distal Offset	L - R .018	L - L .018	L - R .022	L - L .022
1. & 2. Molar - bondable	-25°	0°	0°	205M-36/47	205M-36/47	255M-36/47	255M-36/47
1. & 2. Molar extended base - bondable	-25°	0°	0°	205M-36/47E	205M-36/47E	255M-36/47E	255M-36/47E

MBT prescription

Self-ligating Buccal Tube	Torque	Ang	Distal Offset	U - R .018	U - L .018	U - R .022	U - L .022
1. & 2. Molar - bondable	-14°	0°	0°	206M-16/27	206M-16/27	266M-16/27	266M-16/27
1. & 2. Molar extended base - bondable	-14°	0°	0°	206M-16/27E	206M-16/27E	266M-16/27E	266M-16/27E
Self-ligating Buccal Tube	Torque	Ang	Distal Offset	L - R .018	L - L .018	L - R .022	L - L .022
1. & 2. Molar - bondable	-20°	0°	0°	206M-36/47	206M-36/47	266M-36/47	266M-36/47
1. & 2. Molar extended base - bondable	-20°	0°	0°	206M-36/47E	206M-36/47E	266M-36/47E	266M-36/47E

