

ELASTOMERICS & ACCESSORIES



Chains

90



Elastic Ligatures

91



Intra-oral Elastics

92



Wax

93



Retainer Cases

93

Elastomerics

Elastic Chain | Elastic Ligatures | Separators | Elastic Thread


Elastomerics



Elastomerics

High quality brightly vibrant fun colors inspired by the colors our kids love.

- Medical Grade polyurethane
- Superior color stability, fade resistant
- High elasticity and excellent rebound
- Stain and moisture resistant material



Elastic Chain




Stain and moisture resistant material has consistent force delivery with minimal decay rates.

INTERESTING FACT
A significant reduction in discoloration has been reported with our new Pearl colored Ties and Chains.

Item #	Continuous		Item #	Short		Item #	Long	
EKOAG	silver		EK1AG	silver		EK2AG	silver	
EKOB	blue		EK1B	blue		EK2B	blue	
EKOF	fuchsia		EK1F	fuchsia		EK2F	fuchsia	
EKOG	grey		EK1G	grey		EK2G	grey	
EK0GE	yellow		EK1GE	yellow		EK2GE	yellow	
EK0GR	green		EK1GR	green		EK2GR	green	
EK0HB	sky blue		EK1HB	sky blue		EK2HB	light blue	
EKOK	clear		EK1K	clear		EK2K	clear	
EKOL	purple		EK1L	purple		EK2L	purple	
EK0LC	lavender		EK1LC	lavender		EK2LC	lavender	
EKOO	orange		EK1O	orange		EK2O	orange	
EKOP	pink		EK1P	pink		EK2P	pink	
EKOR	red		EK1R	red		EK2R	red	
EKOS	black		EK1S	black		EK2S	black	
EKOT	turquoise		EK1T	turquoise		EK2T	turquoise	
EKOW	white		EK1W	white		EK2W	white	
EK0MB	metallic blue		EK1MB	metallic blue		EK2MB	metallic blue	
EK0ML	metallic purple		EK1ML	metallic purple		EK0ML	metallic purple	
EKOPearl	pearl		EKOPearl	pearl		EK2Pearl	pearl	

Elastic Ligatures Regular / Glow in the Dark



Designed to maximize treatment efficiency, long-term elasticity facilitate both archwire retention and tooth movement.

Item #	Color		Item #	Color		Item #	color	
ELP-FS	mixed colors		ELP-DG	dark green		ELP-ML	metallic purple	
ELP-FS/F	mixed spring colors		ELP-HB	sky blue		ELP-OGR	metallic ocean green	
ELP-AG	silver		ELP-HR	light pink		ELP-MS	metallic black	
ELP-G	dark grey		ELP-KB	royal blue		ELP-MR	metallic red	
ELP-HG	light grey		ELP-L	purple		ELPG-SF	mixed glow colors	
ELP-US-G	grey ultra-slide		ELP-LL	lavender		ELPG-B	glow blue	
ELP-K	clear		ELP-MP	fuchsia		ELPG-GL	glow	
ELP-US-K	clear ultra-slide		ELP-O	orange		ELPG-GR	glow green	
ELP-S	black		ELP-R	red		ELPG-L	glow purple	
ELP-W	white		ELP-SG	teal		ELPG-O	glow orange	
ELP-AGR	apple green		ELPPEARL	pearl		ELPG-R	glow red	
ELP-DB	navy blue					ELP-MB	metallic blue	
						ELP-MBP	metallic barbie pink	

NOTE: Color range varies from country to country
SIZE: Inner diameter 0.12" 45 sticks per bag

Elastic Thread




Solid elastomeric material extruded in a precision thread form. Exerts continuous long lasting, predictable force. Moisture, odor and stain resistant. Ideal for rotations, retraction and figure 8's.

Elastic Thread order info

Item #	Solid	Item #	Hollow	Item #	Rough
ELF-25G	solid – 025 grey	ELFH-25G	hollow – 025 – grey	ELFU-25G	rough – 025 – grey
ELF-25K	solid – 025 clear	ELFH-25K	hollow – 025 – clear	ELFU-25K	rough – 025 – clear
ELF-30G	solid – 030 grey	ELFH-30G	hollow – 030 – grey	ELFU-30G	rough – 030 – grey
ELF-30K	solid – 030 clear	ELFH-30K	hollow – 030 – clear	ELFU-30K	rough – 030 – clear

Radiopaque Separators




Ultra smooth round edges which makes them easier to insert between the teeth. Our Separators are absorbent free and maintain their elasticity without deterioration.

Radiopaque Separators order info

Item #	Article Description
SSR-2	Separators – radiopaque – blue 100 sticks

Intra Oral Elastics

Latex Elastics



- Force in ounces is measured at 3x passive activation diameter
- Packaged in color-coded tamper-proof patient bags

Diameter	Light 2.5 oz	Medium 3.5 oz	Heavy 4.5 oz	Extra Heavy 6 oz
1/8" / 3.2mm	Red A EG731	Green A EG733	Blue A EG735	Brown A EG736
3/16" / 4.8mm	Red B EG741	Green B EG743	Blue B EG745	Brown B EG746
1/4" / 6.44mm	Red C EG751	Green C EG753	Blue C EG755	Brown C EG756
5/16" / 7.9mm	Red D EG761	Green D EG763	Blue D EG765	Brown D EG766
3/8" / 9.5mm	Red E EG771	Green E EG773	Blue E EG775	Brown E EG776

100 elastics per patient bag: 50 bags per dispenser box - minimum order multiple of 2 boxes

Latex Free Elastics

- Eliminates latex allergic reaction concerns
- Use with aligner cases or ceramic and composite brackets
- Force in ounces is measured at 3x passive diameter
- Packaged in color-coded tamper-proof patient bags

Diameter	Light 2.5 oz	Medium 3.5 oz	Heavy 4.5 oz	Extra Heavy 6 oz
1/8" / 3.2mm	Red A	Green A UNL733	Blue A UNL735	Brown A UNL736
3/16" / 4.8mm	Red B	Green B UNL743	Blue B UNL745	Brown B UNL746
1/4" / 6.44mm	Red C	Green C UNL753	Blue C UNL755	Brown C UNL756
5/16" / 7.9mm	Red D	Green D UNL763	Blue D UNL765	Brown D UNL766
3/8" / 9.5mm	Red E UNL771	Green E UNL773	Blue E UNL775	Brown E UNL776

100 elastics per patient bag: 50 bags per dispenser box - minimum order multiple of 2 boxes

Neon Latex Free Elastics

- Eliminates latex allergic reaction concerns
- Six bright colors in each bag
- Use with aligner cases or ceramic and composite brackets
- Force in ounces is measured at 3x passive diameter
- Packaged in color-coded tamper-proof patient bags

Diameter	Light 2.5 oz	Medium 3.5 oz	Heavy 4.5 oz	Extra Heavy 6 oz
1/8" / 3.2mm	Red A EGF930	Green A EGF932	Blue A EGF934	Brown A EGF936
3/16" / 4.8mm	Red B EGF900	Green B EGF902	Blue B EGF904	Brown B EGF906
1/4" / 6.44mm	Red C EGF910	Green C EGF912	Blue C EGF914	Brown C EGF916
5/16" / 7.9mm	Red D EGF920	Green D EGF922	Blue D EGF924	Brown D EGF926
3/8" / 9.5mm	Red E EGF940	Green E EGF942	Blue E EGF944	Brown E EGF946

100 elastics per patient bag: 50 bags per dispenser box - minimum order multiple of 2 boxes

Patient Wax



- Ideal for fast, temporary relief from fixed brackets and appliances.

Item #	Article Description
DP-WACHS-MULTI	Multi Color / box / 5 strips per box / 50 boxes per pack
DP-WACHS	Clear - box / 3 strips per box / 100 boxes per pack

Retainer Cases



One-piece cases are designed with a locking mechanism that stands up to heavy use and prevents accidental opening resulting in lost retainers.



Item #	Article Description
KD-MB	metallic-blue
KD-MS	metallic-black
KD-MR	metallic-red
KD-MINT	light green
KD-BLAU	box navy blue
KD-GE	yellow
KD-LILA	purple
KD-P	pink
KD-SCHW	black



i ALL METALLIC BOXES CAN BE USED FOR FUNCTIONAL APPLIANCES!