

Tooth	Torque	Tip	Item # .018	Item # .022
11 - Central	17°	4°	106-11	166-11
12 - Lateral	10°	8°	106-12	166-12
13 - Cuspid	-7°	8°	106-13	166-13
13 - Cuspid w hook	-7°	8°	106-13/H	166-13/H
14 - 1. Bicuspid	-7°	0°	106-14/15	166-14/15
14 - 1. Bicuspid w hook	-7°	0°	106-14/15/H	166-14/15/H
15 - 2. Bicuspid	-7°	0°	106-14/15	166-14/15
15 - 2. Bicuspid w hook	-7°	0°	106-14/15/H	166-14/15/H

Item # .022	Item # .018	Tip	Torque	Tooth
166-21	106-21	4°	17°	Central - 21
166-22	106-22	8°	10°	Lateral - 22
166-23	106-23	8°	-7°	Cuspid - 23
166-23/H	106-23/H	8°	-7°	Cuspid w hook - 23
166-24/25	106-24/25	0°	-7°	1. Bicuspid - 24
166-24/25/H	106-24/25/H	0°	-7°	1. Bicuspid w hook - 24
166-24/25	106-24/25	0°	-7°	2. Bicuspid - 25
166-24/25/H	106-24/25/H	0°	-7°	2. Bicuspid w hook - 25

Tooth	Torque	Tip	Item # .018	Item # .022
41 - Anterior	-6°	0°	106-41/42	166-41/42
42 - Anterior	-6°	0°	106-41/42	166-41/42
43 - Cuspid	-6°	3°	106-43	166-43
43 - Cuspid w hook	-6°	3°	106-43/H	166-43/H
44 - 1. Bicuspid	-12°	2°	106-44	166-44
44 - 1. Bicuspid w hook	-12°	2°	106-44/H	166-44/H
45 - 2. Bicuspid	-17°	2°	106-45	166-45
45 - 2. Bicuspid w hook	-17°	2°	106-45/H	166-45/H

Item # .022	Item # .018	Tip	Torque	Tooth
166-31/32	106-31/32	0°	-6°	Anterior - 31
166-31/32	106-31/32	0°	-6°	Anterior - 32
166-33	106-33	3°	-6°	Cuspid - 33
166-33/H	106-33/H	3°	-6°	Cuspid w hook - 33
166-34	106-34	2°	-12°	1. Bicuspid - 34
166-34/H	106-34/H	2°	-12°	1. Bicuspid w hook - 34
166-35	106-35	2°	-17°	2. Bicuspid - 35
166-35/H	106-35/H	2°	-17°	2. Bicuspid w hook - 35