



**Ortho & Aligner**

**Restorations**

**Prophylaxis**

## Indications for controlled minimal invasive enamel reduction treatment

- Slight to moderate crowding
- Bolton Index discrepancies
- Decrease of excessive mesio-inclinations in Upper Class II cases Maintaining lower intercanine dimensions to prevent relapse
- Modification of interproximal contacts to maintain long-term orthodontic stability
- Relapse checking
- Eliminating unaesthetic black triangles in gingival recession cases
- Aesthetic remodelling

- ☐ Order directly to SDC.
- ☐ By your representative in your country. (Ask us for dealer address information)

Article	Description	Item No.	Quantity
<b>SDC IRS ProLign Mini Kit Ortho</b>	<b>Starter Kit ProLign Mini starter-Kit ProLign</b>	<b>1300/Kit 1000/Ortho</b>	_____
<b>Set G5-ProLign</b>	complete ProLign Strip line	<b>1300/S7</b>	_____
<b>Set G5-ProLign</b>	0.10, 0.20, 0.30, 0.40 mm	<b>1300/S4</b>	_____
<b>Refill Packs of 3.</b>	<b>2 side coating</b>		
G5-ProLign	ultra fine 0.10 mm	<b>1310/3</b>	_____
G5-ProLign	fine 0.20 mm	<b>1320/3</b>	_____
G5-ProLign	medium fine 0.30 mm	<b>1330/3</b>	_____
G5-ProLign	medium 0.40 mm	<b>1340/3</b>	_____
<b>Refill Packs of 3.</b>	<b>1 side coating</b>		
G5-ProLign	ultra fine 0.10 mm	<b>1310/1/3</b>	_____
G5-ProLign	fine 0.20 mm	<b>1320/1/3</b>	_____
<b>Refill Packs of 3.</b>	<b>2 side coating</b>		
G5-UltraSoft	ultra fine 0.09 mm	<b>1106/3</b>	_____
<b>Refill Packs of 3.</b>	<b>2 side coating</b>		
G5-CombiStrip	ultra-fine 0.09 mm	<b>1800/UF</b>	_____
G5-CombiStrip	fine 0.12 mm	<b>1800/F</b>	_____
<b>Refill Shank</b>	<b>Universal shank</b>		
G5 Univ. Shank	40 pcs.	<b>1990</b>	_____

Dental office:

Street: \_\_\_\_\_ Tel.: \_\_\_\_\_

ZIP.: \_\_\_\_\_ Fax: \_\_\_\_\_

City: \_\_\_\_\_ Mail \_\_\_\_\_

Dealer stamp



**SDC Switzerland SA**  
Via alla Roggia 2  
CH-6916 Lugano-Grancia  
Switzerland

**T** +41 91 604 54 13  
**F** +41 91 604 54 14  
info@swissdentacare.ch  
www.swissdentacare.ch



Some of the products and designations mentioned in the text are patented or copyrighted. The absence of a special reference or the ® sign should not be interpreted as the absence of copyright protection. This publication is copyrighted. All rights, including those of translation, reprinting, and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm, or other methods) without the editor's prior consent in writing.

Colours and products subject to alterations. Printing errors excepted.



**SWISS**  
**denta**care®  
QUALITY DENTAL PRODUCTS

**SDC Switzerland SA**  
Via alla Roggia 2 • CH-6916 Lugano-Grancia • Switzerland  
T +41 91 604 54 13 • F +41 91 604 54 14 •  
info@swissdentacare.ch • www.swissdentacare.ch

[www.swissdentacare.ch](http://www.swissdentacare.ch)

# Innovative G5 Concept

With the new G5, it is possible to easily exchange the diamond coating by replacing the file and leaving its shank in the head. Attachment and removal of the bend file are made really quick and easy due to the click-mechanism. Once the treatment is completed, the shank can be disposed or, if preferred sterilized to be reused.

The files can be inserted on both sides; by doing so – in case of the diamond-coated files on one side – left as well as right are used with the same file. It is clinically proved that using the G5 Strips in combination with a hand piece, you take 5 time less pression on the teeth.



Reference: 1300/Kit



Reference: 1300/S4



Reference: 1300/S7

## New SDC interproximal calibrated reduction system

### Twin-Aligner Systems and fixed orthodontics

Safe and comfortable treatment. Excellent tools with perfect calibrated thickness, for controlled interproximal reduction. G5-ProLign strips are especially suitable for safe minimal invasive orthodontic treatment.

#### Description

**4 files** precise calibrated thickness (0,10 mm, 0,20 mm, 0,30 mm, 0,40 mm), two-side diamond coated along the entire head length.  
**2 file** one-side-diamond-coating.

#### Indications

- Precise and safe reduction and finishing of proximal enamel surfaces
- High polishing and shaping of proximal enamel surfaces

#### Benefits

- Easier access to the interproximal space with calibrated thickness of the file
- Final high polishing with new extra thin flexible file and diamond coating.
- No iatrogenic damage to adjacent teeth

<b>Name:</b> Orthofile-G5-ProLign	<b>Name:</b> Orthofile-G5-ProLign	<b>Name:</b> Orthofile-G5-ProLign	<b>Name:</b> Orthofile-G5-ProLign	<b>Name:</b> Orthofile-G5-ProLign	<b>Name:</b> Orthofile-G5-ProLign
<b>Thickness:</b> 0.10 mm	<b>Thickness:</b> 0.20 mm	<b>Thickness:</b> 0.30 mm	<b>Thickness:</b> 0.40 mm	<b>Thickness:</b> 0.10 mm	<b>Thickness:</b> 0.20 mm
<b>Grid:</b> ultra fine	<b>Grid:</b> fine	<b>Grid:</b> medium fine	<b>Grid:</b> medium	<b>Grid:</b> fine	<b>Grid:</b> medium fine
<b>Article Ref. 2 sides:</b> 1310/3	<b>Article Ref. 2 sides:</b> 1320/3	<b>Article Ref. 2 sides:</b> 1330/3	<b>Article Ref. 2 sides:</b> 1340/3	<b>Article Ref. 2 sides:</b> 1310/1/3	<b>Article Ref. 2 sides:</b> 1320/1/3
<b>Available packaging refills:</b> 3 pcs.	<b>Available packaging refills:</b> 3 pcs.	<b>Available packaging refills:</b> 3 pcs.	<b>Available packaging refills:</b> 3 pcs.	<b>Available packaging refills:</b> 3 pcs.	<b>Available packaging refills:</b> 3 pcs.
<b>Sterilizable:</b> ▶	<b>Sterilizable:</b> ▶	<b>Sterilizable:</b> ▶	<b>Sterilizable:</b> ▶	<b>Sterilizable:</b> ▶	<b>Sterilizable:</b> ▶
<b>Application:</b> ▶	<b>Application:</b> ▶	<b>Application:</b> ▶	<b>Application:</b> ▶	<b>Application:</b> ▶	<b>Application:</b> ▶



## G5-UltraSoft ultra-thin metal strip with ultra-fine diamond coating

Excellent tool for high polishing or finishing of interproximal areas and subgingival surfaces, will effectively gloss the enamel surface.



<b>Name:</b> Orthofile-G5-UltraSoft	<b>Name:</b> Orthofile-G5-UltraSoft
<b>Thickness:</b> 0.10 mm	<b>Thickness:</b> 0.10 mm
<b>Grid:</b> ultra fine	<b>Grid:</b> ultra fine
<b>Article Ref. 2 sides:</b> 1106/3	<b>Article Ref. 2 sides:</b> 1106/3
<b>Available packaging refills:</b> 3 pcs.	<b>Available packaging refills:</b> 3 pcs.
<b>Sterilizable:</b> ▶	<b>Sterilizable:</b> ▶
<b>Application:</b> ▶	<b>Application:</b> ▶



Patent pending.



<b>Name:</b> G5 Universal Shank	<b>Name:</b> G5 Universal Shank
<b>Compatibility:</b> SIRONA, W&H, NSK, KAVO, SDC	<b>Compatibility:</b> SIRONA, W&H, NSK, KAVO, SDC
<b>Article Ref. 2 sides:</b> 1990	<b>Article Ref. 2 sides:</b> 1990
<b>Available packaging:</b> 40 pieces	<b>Available packaging:</b> 40 pieces
<b>Sterilizable:</b> ▶	<b>Sterilizable:</b> ▶
<b>Application:</b> ▶	<b>Application:</b> ▶